BI 121, Lab 6 Pulmonary Function Testing (PFT)



I. <u>Attendance</u>

II. Pulmonary Function Test/PFT

- A. What? Measure of static & dynamic lung function
- B. Why? Picture of lung health; absence, presence, progression of disease (eg asthma, emphysema); effectiveness of drugs
- C. How? Complete PFT with computer or dinosaur spirometer
- III.<u>Crucial Clinical Measures</u>
 - A. VC vital capacity = FVC forced vital capacity amount of air exhaled after maximal inhalation
 - B. $FEV_{1.0}$ = How much of VC in 1 second? [FEV_{1.0}/FVC] x 100 If \ge 75-80% (0.75-0.80) -> clinically normal
 - If \leq 40-50% –>obstructive disease (eg, asthma)

IV.Your Goals

- A. Estimate your VC & FEV_{1.0} from Nomogram pp 6-6 or 6-7 LLM
- B. Measure these values accurately w/computer PFT LabChart
- C. Compare <u>estimated</u> with actually <u>assessed</u> values to determine whether you're within a healthy range.

Respirometer → measures complete <u>Pulmonary Function</u> <u>Test or PFT</u>!

<u>NB</u>: Should be able to blow out \geq 75 - 85% of VC/FVC in 1 second! That's FEV_{1.0}/FVC \geq 0.75 - 0.85. If less, may indicate asthma or other lung disease.





Normal = Steep



Abnormal = Flatter Downslope (eg, Asthma) *PFT* → measures all lung volumes & capacities (sum of \geq 2 volumes). Subject relaxes & breathes normally into and out of tank.



Sample PFT from Collins 13.5 L Respirometer







1 Estimate, 2 Setup, 3 Assess, 4 Compare

